04-Nov-2014 12:12:17

```
Torsion Suspension (suspTorsion.m)
 - torsion wire diameter (Tungsten, single wire): 466.0401 um - torsion suspension wire length: 0.6 m
 torsion spring constant (2 wire): 0.028786 Nm/radtorsion bar inertia: 0.62492kg*m^2
 - torsion resonance: 0.034158 Hz
You are not injecting squeezing..loozer!
- Seismic Isolator: MinusK
 - Seismic Ground Motion: LLO
 Seismic Isolator: MinusKSeismic Ground Motion: LLO
                               0.100 Watt
Laser Power:
SRM Detuning:
                                0.00 degree
SRM transmission:
                                1.0000
ITM transmission:
                                0.0200
                                1.0000
PRM transmission:
Finesse:
                              313.85
Power Recycling Factor:
                                1.00
Arm power:
                                0.01 kW
                                0.10 W
Power on beam splitter:
Thermal load on ITM: Thermal load on BS:
                                0.000 W
                                0.000 W
Regired TCS efficiency:
                                1.000(estimate, see IFOModel.m for definition)
BNS Inspiral Range:
                               0.000 Mpc
BBH Inspiral Range:
                               0.002 Mpc
Stochastic Omega: 5e+01 Universes
New Nebulous Range:
                               1.093 Mpc
Bar length and diameter: 0.6 m x 0.06 m. (dumbell)
Bar mass (Aluminium, Q~1.13M): 13.911 kg
Bar and Suspension temperature: 293 K and 293 K.
```

